
Infrared microspectroscopy station




The working distance of our microscope is long, and the various attachments can be installed.

•x-y mapping stage minimum step 1 [m
•flow-type cryostat
 Oxford microstat-He
 Temperature 4.2 ~ 400 K
 Window: Quartz, KRS-5, BaF2, polyethylene film
•high temperature DAC
 RT ~ 1000 K, ~30GPa
•low temperature DAC
 10 ~400 K, ~20GPa
 Window: Quartz, KRS-5, BaF2, polyethylene film

Infrared microspectroscopy station Far infrared region

Air transmission spectra

Si-bolometer Detector Resolution 2 cm⁻¹ Accumulation number 0.6 0.5 0.4 900 0.3 0.9 0.0 400 600 Wavenumber (cm⁻¹) Wavenumber (cm⁻¹) SR = gl x 20 @300 cm⁻¹ noise @300 cm-1 SR:5% gl:10%




Infrared microspectroscopy station Mid infrared region

 Air transmission spectra


 Detector
 MCT

 Resolution
 4 cm-1


 Accumulation number
 124



Beam profile at the focal point of the microscope through 2 m pinhole (700 ~9000 cm-1)



Results



Absorption spectra of alumina polycrystals

Heat treatment of alumina hydrate produces a series of transition aluminas, such as []-, []-, and so on, before final conversion to the most stable []-phase Al2O3.

(a) □-phase, ~30 □m diameter (b) □-phase, ~ 40 □m diameter

Magneto-optical spectroscopy station

Spectral range	800 ∼ 20000 cm ⁻¹
Magnetic filed	~14 T
Temperature	~4 K
Spatial resolution	10 □m

